



## MICRO MOSQUITO FORCEPS

Used as hemostats to control bleeding in  
blood vessels.

## MICRO MOSQUITO FORCEPS

Our micro mosquito forceps have a very slender pattern with vertical serrations on the jaws and are typically used as hemostats to control bleeding in blood vessels. The slender pattern of these forceps indicate that they are lightweight and allow the doctor to possess control of the instrument.



- DESIGNED FOR CLAMPING SMALL DELICATE BLOOD VESSELS AND TISSUE
- RATCHETED LOCKING MECHANISM
- CONTROLS HEMOSTAT'S JAW PRESSURE
- DELICATE PATTERN
- FULLY SERRATED JAWS

**SLENDER PATTERN WITH VERTICAL SERRATIONS.  
IT CAN HOLD THE BLOOD VESSELS AND PREVENT BLOOD LOSS  
IN SURGICAL PROCEDURES.**





**ASPIRE**  
GLOBAL HEALTHCARE

[www.aspireintl.co](http://www.aspireintl.co)

Office 2108, Bayswater Tower - Marasi Dr - Business Bay - Dubai.